

## Film, CR Cassette, and DR Cassette Protection

**No more broken or damaged cassettes!** The unique design of the Cassette Protector provides an easy way to take weight-bearing x-rays of feet (AP view). The Cassette Protector is simply placed over a cassette which is first positioned on the floor. Patients step onto the Cassette Protector with as much weight as needed to get the desired image.

- Optimal positioning for weight-bearing AP foot images, of even the heaviest patients
- The face plate is made of unbreakable polycarbonate with designs supporting 275 – 400 pounds
- Keeps patient's foot as close to the target image platform as possible to ensure high-quality images with no significant magnification
- The face plate is X-ray lucent, allowing the X-rays to pass through the Cassette Protector with very little absorption or scattering
- The non-slip rubber floor grips keep the cassette protector from slipping on a hard floor
- Easy to keep clean, the Cassette Protector is made from germ-resistant nonporous polyethylene, polycarbonate, and stainless steel
- Simple, reliable, and virtually indestructible
- Cassette Protectors are available for all sizes and models of film, CR, and DR imaging systems



10 x 12 inch or 24 x 30cm

**Item 2934: 275 lb Capacity - \$399**  
**Item 2934H: 450 lb Capacity - \$474**



14 x 17 inch

**Item 3947: 275lb Capacity - \$489**  
**Item 3947H: 450lb Capacity - \$589**

*Clear Image Devices' series of positioning devices is the result of innovative research and design, responding to the expressed needs of radiographers. CID produces only high-quality products that enhance medical care and promote the continued advancement of medical imaging.*